English

SUPLASYN® SUPLASYN® 1-SHOT

STERILE SODIUM HYALURONATE SOLUTION

HYALURONIC ACID (HA) is a normal component of the synovia and plays a central role in maintaining the physiological internal environment of the joint.

COMPOSITION: Viscoclastic solution of a defined molecular weight of purified hyaluronic acid, produced by fermentation.

Each syringe contains:

	SUPLASYN	SUPLASYN 1-SHOT
Hyaluronic acid sodium salt	20 mg	60 mg
Excipients q.s	2 ml	6 ml

PROPERTIES: Hyaluronic acid is the prototype of a wide range of saccharide biopolymers (glycosaminoglycans), important components of all extracellular tissue structures, including cartilage and synovial fluid. The active substance of SUPLASYN is a hyaluronic acid of defined molecular chain length with a high degree of purity. The introduction of SUPLASYN into the synovial space will assist in the normalization of the joint following arthrocentesis.

INDICATIONS: SUPLASYN and SUPLASYN 1-SHOT are indicated for the symptomatic treatment of osteoarthritis. SUPLASYN has been shown to be beneficial in osteoarthritis for the management of pain and improvement in physical function of joints. More than one joint may be treated at the same time.

DOSE AND ADMINISTRATION: Depending upon joint size, up to 6 ml may be administered intra-articularly.

SUPLASYN: The recommended schedule for SUPLASYN supplied in a 2ml syringe is 1 injection per week for 3 weeks, but up to 6 may be given depending on patient's condition.

SUPLASYN 1-SHOT: Intended for single administration.

Use strict aseptic technique. Discard any unused portion of the syringe. To use the pre filled syringe, remove the Luer lock cap, attach a suitable cannula (recommended is 21 – 25 G depending on joint) and secure it by turning slightly. GRADUATION ON THE SYRINGE LABEL IS TO BE USED AS A GUIDE ONLY.

CONTRAINDICATIONS/PRECAUTIONS: Do not administer to patients with known hypersensitivity reactions. Respect usual precautions and contraindications for any intra-articular injection. Do not inject intra-vascularly. SUPLASYN should not be used in patients presenting an inflammation/irritation of the joint, since adverse events more commonly occur in patients with already existing joint inflammation/irritation. As no clinical evidence is available on the use of Hyaluronic Acid in children, pregnant and lactating women, treatment with SUPLASYN is not recommended in these patients. The patient should rest 24-48 hours after the injection and avoid any strenuous activity over the full course of the treatment.

Transient short duration pain may occur following intraarticular introduction. The affected joint may show a mild local reaction like pain, feeling of heat, hyperthermia, redness, effusion, irritation, and swelling/inflammation. If these symptoms occur, rest the affected joint and apply ice locally. Symptoms subside within days for most of the patients. In some cases, mild local reactions such as pain, irritation, swelling/joint inflammation and effusion may be significantly enhanced and much more severe as an expression of hypersensitivity. In such cases, a therapeutic intervention could be necessary, e.g. aspiration of joint fluid. Local adverse reactions could be accompanied by systemic reactions such as fever, chills, or cardiovascular reactions, and in rare cases anaphylactic reactions. In extremely rare circumstances, rash/itching, urticaria, synovitis, and a drop in blood pressure have been reported following the administration of SUPLASYN. Discontinue use if adverse reactions are experienced. Avoid using SUPLASYN with sterilising or sanitising agents containing quaternary ammonium salts solutions.

WARNING: KEEP OUT OF REACH OF CHILDREN. DO NOT USE IF BLISTER IS DAMAGED. TO BE USED BY A PHYSICIAN ONLY. FOR SINGLE USE ONLY. DISCARD UNUSED PORTION OF SYRINGE. REUSE MAY CREATE A RISK OF CONTAMINATION AND/OR CAUSE PATIENT INFECTION OR CROSS INFECTION.

PACKAGING: Available in 2ml or 6ml syringes; packages of 1 and 3.

STORAGE: Store between 4°C and 25°C. DO NOT FREEZE.

Bring to room temperature before injection.

Date of preparation: February 2017

Date of revision: First Issue

Mylan Institutional, Inverin, Co. Galway, Ireland.
© Mylan Pharma Group Limited, Ireland

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Mylan Institutional, Inverin, Co. Galway, Ireland.

© Mylan Pharma Group Limited, Ireland

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
Operator **Yarnold, Paul**

Date **20/02/2017 15:53:46**

Attention: The colours displayed in this proof are for electronic representational purposes only. Please refer to the PANTONE® colour guides for accurate colour matching. By approving this file, you are signifying full approval of design and text. To ensure accurate viewing and printing of this PDF file, Overprint Preview must be activated for accurate colour representation. MPS Leaflets PDF files are designed to work with Adobe Acrobat 6 or higher and are PDF/X -01a-2001 compliant.

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Operator Yamold, Paul

Date 20/02/2017 15:53:46

Attention: The colours displayed in this proof are for electronic representational purposes only. Please refer to the PANTONE® colour guides for accurate colour matching. By approving this file, you are signifying full approval of design and text. To ensure accurate viewing and printing of this PDF file, Overprint Preview must be activated for accurate colour representation. MPS Leaflets PDF files are designed to work with Adobe Acrobat 6 or higher and are PDF/X -1a-2001 compliant.